



AMX 4 portable x-ray

# GE AMX 4 Portable

## Remanufactured AMX 4



GE AMX 4 Portable

**Huestis Medical remanufactures all GE AMX portables under a thorough electronic, mechanical and cosmetic program. Numerous upgrades include a new maintenance-free battery pack, critical structural reinforcements and replacement of major cabling and control components. Huestis Medical then calibrates and tests each system to exacting standards backed by an 12-month warranty on workmanship and parts.**

**Ask about our extensive inventory of new replacement battery systems and new, used and refurbished AMX parts.**

### Features

- New 28 Amp/Hour Lead Calcium Batteries.
- New control panel decals.
- New H.S. cord and housing.
- New small diameter HV cables.
- New L.V. cables.
- New counterweight cables.
- New wrap around bumpers.
- Upgraded front casters (if needed).
- New drive belts.
- New hospital grade cord cap.
- New collimator dials.
- Mechanical assemblies inspected and rebuilt as necessary.
- New bearings where needed.
- Complete calibration to meet or exceed OEM specifications. [AMX-4 UPGRADES:]
- Reinforced Base
- New style spring loaded casters
- New style cassette draw
- Pad added to the bottom of the cassette drawer
- Hall effects tube latch assembly
- New style relay assembly for tube latch
- Upgraded drive wheels
- Upgraded drive motors
- New style micro-switches in the drive handle
- New style micro-switches in the in the collimator handles
- New style knobs and pointers on collimator
- Proprietary upgrade to rotor board
- Proprietary upgrade to the CPU
- ESD (Electro Static Discharge) cable added to the base unit
- ESD wire added to the top shroud
- ESD added to the drive handle
- New style PCB interconnect cables with choke installed

- Protective wrap added to the top shroud
- AMX-4+ style vertical column positioning of xray tube through 270-degree column rotation
- Proprietary battery charge labeling
- Green charge indicator light added
- Remote battery charge voltage test terminal block added
- A copy of the calibration data (battery back-up RAM chips) is archived for a period of 2 years

### Technical Specifications

- 100 mA constant.
- kVp 50-125.
- mAs 0.4-320.
- 24 kVp stations.
- 30 mAs stations.
- X-Ray tube MX-75 275k HU.
- Focal spot .75mm.
- Dual motor drive (forward & reverse).
- Variable speed drive.
- High frequency generator.
- LED display.
- Electric solenoid locks.



### AMX-4 - High Frequency Performance

The GE AMX-4 features a high frequency, microprocessor controlled generator and are upgraded to a vertical column. These thoroughly remanufactured units are backed by a 12 month warranty.



### Replacement Battery Upgrades

Our sealed, Lead-Calcium battery packs offer maintenance-free service without modification. Ask about our retro-fit kits.

**HUESTISMEDICAL**  
making it possible™

Huestis Medical, 597 Winthrop Street  
Taunton, Massachusetts 02870-2166 USA  
508-824-8628 • 800 333-6070  
Fax 508-824-8629  
www.huestismedical.com

ISO9001:2000  
REGISTERED