

## GE AMX 4 Digital Remanufactured Portable With Canon Digital Upgrade

Upgrade your existing AMX 4 or purchase an AMX 4 with the digital upgrade already installed. You can have a mobile DR at your fingertips by retrofitting to either a CXDI-50G or a CXDI-31 panel. The speed of the Canon Digital Radiographic system allows you to see your study right from the portable, brightly displayed on the easy to-use touch screen.

### AMX 4 with Canon Digital Upgrade Specifications

#### CXDI-50G Digital Receptor

- Manufacturer - Canon
- Model Name - CXDI-50G
- Detector Array - Flat Panel, a-Si
- Active Image Size - 14 in x 17 in (35 cm x 43 cm)

#### CXDI-31 Digital Receptor

- Manufacturer - Canon
- Model Name - CXDI-31
- Detector Array - Flat Panel, a-Si
- Active Image Size - 9 in x 11 in (22.6 cm x 28.8 cm)



GE AMX 4 Portable  
with digital panel upgrade



**ISO9001**  
REGISTERED