



GE RFX Package Like New RFX Suite with New X-ray Generator

Huestis Medical's like new GE RFX systems begin with a complete tear-down. All mechanical and electronic systems are then functionally and cosmetically inspected, evaluated and remanufactured to exacting standards. After a detailed three-step paint process, each system is pre-staged for operational compatibility and calibration to OEM specifications.

Bundled with a NEW high frequency x-ray generator from CPI and an array of powerful imaging options, this versatile suite offers excellent imaging at an affordable price.

Features:

Table Features

- Quiet angulation drive with soft start and stop feature.
- Tabletop longitudinal drive interlocked with angulation drive to prevent collision with floor and ceiling.
- Tabletop and footrest extend downward toward floor with table vertical.

Range of Motion

- 73.5 in (187 cm) fluoroscopic beam center-to-floor dimensions when table is vertical for cervical esophagus coverage on patients up to 6 ft 8 in (2.0 m) tall.
- 11 in (28 cm) transverse fluoroscopic carriage travel.
- 20 in (50 cm) maximum caliper opening between bottom of spotfilmer or pedestal and tabletop, 19 in (48 cm) maximum with 4-way flat top.
- Structural plastic front and end panels shaped to eliminate sharp edges and barium traps while resisting strains, dents and corrosion.
- Undertable x-ray tube equipped with a forced oil/air heat exchanger for improved cooling.

Convenient Operator Features

- All table controls needed by technologist are located conveniently on table.
- Interlocked patient step eliminates need for accessory foot stool.
- Patient hand grips may be quickly attached at any point on front or rear of tabletop.

Myelographic Stop

- When automatic stop at horizontal is defeated for myelography and other needle intensive procedures, the compression cone drive is disabled.

Efficient Bucky Design

- Entrance type, three field Quantamat Ionization Bucky Detector designed for minimum object film distance.
- Grid lines virtually invisible at 8 m/sec with 103-line grid.

RFX 88 Spotfilm Device

- Accepts 9.5 in x 9.5 in (24 x 24 cm) cassette sizes for operating simplicity and rapid access from fluoro mode.
- Effort sensing, longitudinal power assist control in positioning handle.
- 1-on-1, 2-on-1 longitudinal, 2-on-1 transverse, and 4-on-1 formats.
- Spotfilm access time from fluoro mode is approximately 0.75 seconds.
- Rapid-sequence spotfilms are accomplished with approximately 0.5 seconds between exposures.
- Table Bucky "center-finder" feature.
- 88/35 or DSP also available.

Included Packaged Components

- Like New GE RFX table and GE RFX 88 Spotfilm Device.
- New CPI high-frequency generator.
- Like New GE XT Suspension and GE Sentry III x-ray Collimator.
- 9/6/4.5 in I.I. (12 in, 14 in 16 in available)



CPI High Frequency Generator

CPI incorporates time tested technology with innovative user interfaces to meet the most demanding imaging needs.



Digital Imaging at Your Fingertips

Optional new digital imaging system, provides user-friendly, crisp fluoroscopic/spot imaging.

Huestis Medical/ARI 68 Buttonwood Street, Bristol, Rhode Island 02809 USA
phone 401 253 5500 800 972 9222 fax 401 253 7350 www.huestis.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA


Best

healthcare for everyone

OPTIONS



Like New GE XT Suspension

GE XT Suspension systems are structurally inspected and thoroughly remanufactured to provide accurate radiographic positioning.



Full Line of Like New R/F Suites

We offer a variety of complete, packaged room systems with versatile system options. Ask about our complete like new line.



Optional Touch Screen Available

CPI high-frequency x-ray generators offer unsurpassed power and reliability. Extremely versatile, they are also a perfect match for digital imaging upgrades.



Like New Sentry III Collimators

GE Sentry III collimators automatically sense the size of the cassette and the source-to-image distance (SID).



Optional Large-Format Image Intensifier

Optional 12 in, 14 in, or 16 in I.I. offers excellent imaging results at an affordable price for many large format applications.

Like New RFX 88 Spotfilm Device

GE RFX 88 spotfilmer has automatic loading of 9.5 in x 9.5 in (24 x 24 cm) cassettes from either the front or rear position.



Thoroughly Remanufactured RFX Table

Our remanufactured R/F systems undergo a thorough cleaning, structural inspection, component replacement, refinishing, pre-staging and calibration to OEM specifications.



Affordable Digital Imaging

Optional Digital R/F imaging packages make exceptional digital imaging available to fit the typical budget.

GE RFX Package Included Components

This comprehensive system includes our most popular like new and new R/F imaging components bundled into a single "turnkey" package. Pre-staged for system compatibility and calibrated to OEM specifications, this versatile suite is a completely viable, and reliable alternative to new equipment at about half the cost. Ask about our complete line of remanufactured imaging systems tailored to your imaging needs.

