

## Huestis Fixed Lasers — Red or Green Diode

- Unique three-filter design reduces transient power issues and extends diode life
- Precise line quality to ensure accurate position reproducibility
- Easy access to mounting bolt during installation
- Stable aluminum housing

	Red Diode Laser	Green Diode Laser
Depth	3.9" (9.9 cm)	
Width	3" (7.6 cm) housing; 5.5" (13.9 cm) mounting plate	
Height	11.5" (29.2 cm) housing; 13" (33 cm) mounting plate	
Laser Beam Output	Less than 1 mW per beam (Class II laser)	
Line Width	0.9–1.2 mm at isocenter when mounted 3 m from isocenter	0.9–1.1 mm at isocenter when mounted 3 m from isocenter
Fan Beam Angle	Greater than 50 degrees	
Wavelength	635–650 nm	520 nm

