

HUESTISMEDICAL

making it affordable™



Huestis Medical/ARI 68 Buttonwood Street, Bristol, RI 02809 USA
phone 401 253 5500 800 972 9222 www.huestis.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA



Best®

healthcare for everyone

BMI.v.08.2010

Huestis Medical

Remanufactured, Not Refurbished

What is the difference?

Remanufactured

New Generators

New Cameras

New Cables

Digital upgrades available

All equipment and components are thoroughly dismantled, inspected, cleaned and repaired or replaced

Cosmetic upgrades

Fully pre-staged and calibrated

Refurbished

Cleaned

Fresh coat of paint

Any other cosmetic upgrades

Huestis Medical/ARI has over 20+ years
experience in remanufacturing

Quality you can count on

ISO 9001 Registered Company



HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

RFX

RFX PACKAGE – RADIOGRAPHIC/FLUORO SYSTEMS

- NEW X-ray Tubes
- NEW Image Intensifier
- NEW Generator
- NEW Safety cables, high & low voltage cables and all new hardware
- Digital upgrades available
- Fully pre-staged and calibrated

**80% NEW
20% ORIGINAL**



NEW Dunlee Tube



HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

RFX

R/F and Rad Units Come
with a Standard 2-year Warranty



Digital Imaging – At Your Fingertips



GE XT Suspension



Touchscreen

DID YOU KNOW?

All Huestis Medical remanufactured systems are guaranteed to deliver performance and safety features that meet or exceed original manufacturer's specifications – at about half the price of new.



88/35 Spotfilm Device

HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

Hydradjust® IV

HYDRADJUST® IV – REMANUFACTURED UROLOGY TABLE



- Smooth and reliable hydraulic functions
- Superior patient positioning and comfort
- Custom matched to your specific needs
- NEW High-frequency generator

Contact us for detailed product brochure



HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

AMX IV

AMX IV PORTABLE – REMANUFACTURED AMX IV

**12 MONTH
WARRANTY**



- NEW Batteries
- NEW Cables
- NEW Housing, belts, knobs and much more
- ▶ 270° Column rotation



Remote Tech Switch

**AMX III, AMX II AND AMX 110
ALSO AVAILABLE!**

**18 MONTH
WARRANTY**



HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

MINI C-ARM

FLUOROSCAN MINI C-ARM – EXTREMITY IMAGING

- ▶ 6/4 dual mode image intensifier
- High intensity monitor
- Smooth 360° rotation
- Smallest available focal spot x-ray tube; Unrivalled resolution and clarity with significant reduction in radiation exposure
- Unique command panel; Puts seven essential commands into the sterile field
- ▶ 6/4 Dual Mode; Selectable field size for superior visualization
- User focused; Designed to be efficient and easy to use



We Will Not Be Undersold!


**12 MONTH
WARRANTY**

HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

OEC 9800

OEC 9800 C-ARM – REMANUFACTURED OEC 9800 C-ARM



**12 MONTH
WARRANTY**

- Restored to like-new condition
- 1k x 1k High resolution CCD technology
- Gold standard and mobile fluoroscopy applications



Dicom options are available for OEC 9800, OEC 9600, and OEC 9400.

HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

OEC 9600

OEC 9600 C-ARM – REMANUFACTURED OEC 9600 C-ARM

**12 MONTH
WARRANTY**

- Offered with new and upgraded imaging components
- Units disassembled, mechanically restored, refinished, calibrated and pre-staged prior to delivery
- High frequency generator
- ESP, vascular and neuro vascular modules available



**Guaranteed
Performance**



HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

OEC 9400

OEC 9400 C-ARM – REMANUFACTURED OEC 9400 C-ARM

- Offered with new and upgraded imaging components
- Units disassembled, mechanically restored, refinished, calibrated and pre-staged prior to delivery
- High frequency generator

**12 MONTH
WARRANTY**



All of our equipment guaranteed to meet or exceed OEM specifications. Contact us for more detailed product brochures.

HUESTISMEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com

PARTS

PARTS

We carry parts for all our remanufactured equipment:

R/F • Rad • Urology • GE AMX • OEC C-Arm



Tubes




Leg Supports



Circuit Boards

We have added to our already extensive inventory with our newest division – Best Medical Parts. We now carry parts for GE, Phillips, Toshiba, Shimadzu, Siemens and many other companies.

Give us a call at 866 602 8846 or 800 972 9222.



Replacement Battery Upgrade

Our sealed, Lead-Calcium battery packs offer maintenance-free service without modification.



AMX Replacement Parts

Our extensive new, used & refurbished AMX parts inventory is an excellent resource for quality replacement parts.

For AMX or C-Arm parts or batteries contact:
800 747 2786 or
847 426 1055

HUESTIS MEDICAL

phone 401 253 5500 800 972 9222 www.huestis.com www.teambest.com